## Rules for storing samples in a walk-in fridge or freezer

	<ul> <li>Consolidate your samples into a storage container</li> <li>Use a cooler, plastic tote or strong cardboard box</li> <li>Match the volume of outer container = contents (don't waste space)</li> <li>Ensure box is strong, bottom is taped and not too large</li> <li>Containers should be able to fit on the shelves</li> <li>All containers must have lids</li> </ul>
Ŕ	<ul> <li>Label each container (top and sides)</li> <li>We should be able to determine who the container belongs to without pulling it off shelf/having to lift it up</li> <li>Label with:         <ul> <li>Your 1st/last name</li> <li>Supervisor's name</li> <li>Date</li> <li>Remove all other names/addresses if reusing container</li> <li>Where possible write directly on container with permanent marker. Tags/tape fall off</li> </ul> </li> </ul>
	Absolutely no storage of lab chemicals of any kind e.g. Sigma, Fisher chemicals. These must be stored in your lab NO storage of flammable samples. Use a spark free fridge/freezer specific for flammable storage
	All animal/fish samples must be stored in bags that are leak proof Ideally in hard sided container NO exposed tissue, fur NO exposed bones, horns –nothing than can scratch/cut
	Remove & discard samples when your research is complete Ensure you have obtained your supervisor's approval Discard appropriately e.g. garbage, biohazardous/hazardous waste
	Unlabelled samples may be discarded If we are unable to ascertain owner If fridge/freezer fails and we have to relocate contents temporarily

If you have any questions, contact Science Facilities:

Angela Sikma angelasikma@trentu.ca or Chris Williams cwilliams@trentu.ca